ÖSSUR
LIFE WITHOUT LIMITATIONS

Navii[®]

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Table 1 - Warning Signals

Warning Type	Device State	Device Feedback	Action
Low battery	Low battery (5% charge)	Beeps 10 times every 7 seconds. Vibration 3 times every 10 seconds. Device status indicator turns red. One battery indicator flashing green.	Charge battery.
	Critical low battery (2% charge)	Automatic shutdown.	Charge battery.
Device temperature	Critical high device temperature	Beeps 10 times every 7 seconds. Vibration 3 times every 10 seconds. Device status indicator turns red. Automatic shutdown after 2 minutes.	Stop using device. Send to Össur for service.
Critical error	Critical hardware error	Beeps 10 times every 7 seconds. Vibration 3 times every 10 seconds. Device status indicator turns red.	Stop using device. Send to Össur for service.
	Forced shutdown error	Automatic shutdown.	Stop using device. Send to Össur for service.

Automatic shutdown
 When the battery is depleted, the device encounters a critical error or high temperature, the device will emit a long steady beep with pulsing vibration for 30 seconds while the status indicator turns red.
 The device will then emit three descending tones, vibrate three times, and shut down.

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Unity Elevated Vacuum

Seal-In Suspension

The benefit of using Seal-In liners compared to other kinds of suspension systems is supported in:

- Secure suspension¹
- Reduced pistoning^{1,2,3,4,5}
- Improved rotational control³

1. Eshraghi, A. et. al. 2014

2. Brunetti, S. et. al. 2013

3. Ali, S. et. al. 2012

Elevated Vacuum Suspension

Adding an elevated vacuum system to a prosthesis with a Seal-In liner leads to:

- Improved transitions in body positions⁵
- Improved gait symmetry^{4,9}
- Improved balance and stability⁶
- Reduced volume fluctuations^{4,10,11}
- Reduced risk of falls^{4,7}

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


Navii® Leg with Unity

For Unity TF Kit used with NAVii® knee and cover:

Note: The Tube should be placed on the medial side of the prosthesis.

1. Stamp two 8mm holes on the medial side of the knee cover (Fig. 6).
Use the small circular indentations on the knee cover as placement guides. (8mm drill)
2. Attach the proximal and distal tube clips to the knee cover by inserting the tube clips through the holes and securing them with the rivets (Fig. 7).
3. Thread the tube through, then place the small tube fixation on the pylon above the knee (Fig. 8a).
4. Place the tube clamp on the pylon between foot and knee (Fig. 8b).
5. Thread the tube through the tube clamp (Fig. 8b) and the two tube clips (Fig. 8c).
6. Make sure the tube w/protective sleeve is placed next to the knee joint and not behind it (Fig. 8d).

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
Lisan Scheepers
Clinical Specialist, Physiotherapist

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

Rehabilitation After Amputation: Paqueras?
 Impairment variables producing activity limitation in individuals with lower limb amputation. Ray
 The influence of falling, fear of falling, and balance confidence on prosthetic mobility and social activity among individuals with a lower extremity amputation. Miles

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Collingham TS, Peckitt LS, Shore AD. Reamputations, Mortality, and Health Care Costs Among Persons with Dysvascular Lower Limb Amputations. *Arch Phys Med Rehabil*. 2005;86:180-6.

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
Rehabilitation After Amputation: Experience
Amputation enables predicting activity limitation in individuals with lower limb amputation, Riva
The influence of falling fear, fear of falling, and balance confidence on prosthetic mobility and social activity among individuals with a lower extremity amputation, Miller

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Percent increase in energy expenditure

Esquenazi and DeGiacomo: "Rehabilitation After Amputation," Journal of American Podiatric Med. Assoc. 61(1): 13-22, 2001

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Rehabilitation After Amputation. Esquenazi
 Impairment variables predicting activity limitation in individuals with lower limb amputation. Rays
 The influence of falling, fear of falling, and balance confidence on prosthetic stability and social activity among individuals with a lower extremity amputation. Miller

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WHO. Global Recommendations on Physical Activity for Health. Geneva: World Health Organization; 2010.

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National Sporting Goods Association Sport Participation in 1990. MLI Prospect, NYSGA, 1991

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1. Health Council of the Netherlands. *Sedentary behaviour and risk of chronic diseases*, 2017
2. J. Y. Chau, et al. Daily Sitting Time and All-Cause Mortality: a Meta-analysis. *PLoS ONE*. 2013; 8(11)

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Rehabilitation After Amputation, Esquerazi

Impairment variables predicting activity limitation in individuals with lower limb amputation, Raya

The influence of falling, fear of falling, and balance confidence on prosthetic mobility and social activity among individuals with a lower extremity amputation, Miller

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Galley, R. (2008). Secondary conditions related to prosthetic users and 10 steps to reduce the risk of injury.

Strauß, Peter A., et al. "The prevalence of osteoarthritis on the intact hip and knee among traumatic leg amputees." Archives of physical medicine and rehabilitation 90.3 (2009): 440-448.

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
Armstrong et al. (2011) Journal of Bone and Mineral Research, Volume 26, Issue 6

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² D. M. Ehde, et al. Chronic phantom sensations, phantom pain, residual limb pain, and other regional pain after lower limb amputation, *Arch Phys Med Rehabil* 2000; 81(10):9-1044.

³ D. M. Ehde, et al. Back pain as a secondary disability in persons with lower limb amputations, *Arch Phys Med Rehabil*, 2001, 82:731-4.

⁴ D. C. Smith, et al. Phantom limb, residual limb, and back pain after lower extremity amputations, *Clin Orthop*, 1999, 361:29-38.


⁵ Kulikarni, et al. Chronic low back pain in traumatic lower limb amputees, *Clinical Rehabil*, 2005, 19(1):81-86.

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
Goujon-Pillet, H. Et al. Three dimensional motions of trunk and pelvis during transfemoral amputee gait. Arch Phys Med Rehabil. 2008;89(1):87-94.

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
1. Deyo RM, Weinstein JD. Low back pain. N Engl J Med. 2001;344:363-9.

2. Junk S, Terjesen T, Rossvoll I, Braten M. Leg length inequality measured by ultrasound and clinical methods. Eur J Radiol. 1992;14:185-8.

3. Defrin R, Benaymin S, Aldubi R, Pick C. Conservative correction of leg-length discrepancies of 10mm or less for the relief of chronic low back pain. Arch Phys Med Rehabil. 2008;89:2675-80.

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Rehabilitation After Amputation: Esauvaz

Important variables predicting activity limitation in individuals with lower limb amputation: Pajo

The influence of feeling, fear of falling, and balance confidence on proximal mobility and locomotion among individuals with a lower extremity amputation. Milg

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Miller WC, Deathe AB, et al. The influence of falling, fear of falling, and balance confidence on prosthetic disability and social activity among individuals with a lower extremity amputation

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Miller WC, Deathe AB, et al. The influence of falling, fear of falling, and balance confidence on prosthetic mobility and social activity among individuals with a lower extremity amputation

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Rehabilitation After Amputation: Experiences
Impairment variables predicting activity limitation of individuals with lower limb amputation: Ayo
The influence of falling, fear of falling, and balance confidence on prosthetic mobility and social activity among individuals with a lower extremity amputation. Miller

Dual-task standing and walking in people with lower limb amputation: Sara J Morgan, Brian J Hafner, Deborah Kartin and Valerie E Kelly Prosthetics and Orthotics International 1-15 Article first published online July 26, 2018

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
Rehabilitation After Amputation: Equilibrium
 Improvement variables predicting activity limitation in individuals with lower limb amputation. *Ally*
 The influence of feeling, fear of falling, and balance confidence on prosthetic mobility and social activity among individuals with a lower extremity amputation, Miller

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- Too Long
- Too heavy
- Pain

Oliver A, Helgmann D, Ruckw D, Barnett C, Milson P, Scheepers L, Berh M & Miller WC. A Meta Analysis of Training Interventions to Improve Gait, Postural Control, Falling, and Fear of Falling in People With a Lower Limb Amputation. Poster presented at the Canadian Association of Physical Medicine and Rehabilitation (CAPMRE) Annual Scientific Meeting, Montebello, QC, Canada.

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