
743A200 PROS.A. Assembly

DE	Gebrauchsleitung (Benutzer / Techniker)	3
EN	Instructions for Use (User / Technician)	21
FR	Mode d'emploi (Utilisateur / Technicien)	39
ES	Instrucciones de uso (usuario / técnico)	57
DA	Brugsvejledning (bruger / tekniker)	75
TR	Kullanım kılavuzu (Kullanıcı / Teknisyen)	93
RU	Руководство по применению (для пользователей / технического персонала)	109

DE

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

PROS.A. Assembly Instructions for Use

Contents	Page
1 General Information	23
1.1 Foreword	23
1.2 Intended Use	23
1.3 Function	23
1.4 Connection Systems	24
1.5 Liability	24
2 Safety Instructions	24
2.1 Explanation of Symbols	24
2.2 General Safety Instructions	24
3 Delivery and Preparation for Use	26
3.1 Instructions for Safe Transportation	26
3.2 Unpacking and Checking	26
3.3 Place of Installation and Environmental Conditions	26
3.4 Installing the Alignment Apparatus	26
4 Alignment of a Modular Transfemoral Prosthesis	27
4.1 Foot Positioning	27
4.2 Knee Joint Positioning	28
4.3 Socket Positioning	30
4.4 Documenting the Prosthesis Alignment	33
4.5 Using a Laser	34
4.6 Static and Dynamic Alignment Optimisation	35
5 Alignment of Knee Disarticulation Prostheses	36
6 Maintenance and Cleaning	37
6.1 Cleaning	37
7 Technical Data	37
8 Disposal Information	37


Component Overview

- 1 Holder for adapter inserts and socket clamping fixtures
- 2 Mounting frame for the bench alignment of leg prostheses
- 3 Socket adapter height support
- 4 Folding crank for the hydraulic adjustment of the socket adapter height support
- 5 Socket clamping fixture
- 6 Knee bracket for knee joint adapter inserts

- 7 Knee joint adapter height support with clamping levers
- 8 Adjustment plate for heel height adjustment
- 9 Footplate with angle adjustment
- 10 Scale for the adjustment of the knee-floor measurement
- 11 Scale for the adjustment of the trochanter-floor measurement
- 12 Height-adjustable feet

- 1 Knee joint adapter inserts
- 2 Socket clamping fixture
- 3 Centring unit extension

- 1 Centring unit

INFORMATION

Date of the last update: 2014-03-12

- Please read this document carefully.
- Follow the safety instructions.

1 General Information

1.1 Foreword

With the acquisition of the PROS.A. Assembly alignment apparatus (**PRO**sthesis **A**lignment **A**ssembly), you have selected a quality product for the bench alignment of modular leg prosthesis systems. Before you start using the alignment apparatus, please make sure you read the section entitled “Delivery and Preparation for Use” as well as the safety information.

Section 4, “Alignment of a Modular Transfemoral Prosthesis”, illustrates how to position the components of a transfemoral prosthesis in the alignment apparatus and explains the options to document the prosthesis alignment. Section 5 describes particularities for the alignment of knee disarticulation prostheses in the PROS.A. Assembly. Please contact Ottobock if you have any questions or concerns that go beyond these instructions for use.

Technical changes to the design described in these instructions for use are reserved.

1.2 Intended Use

The Ottobock PROS.A. Assembly alignment apparatus is suitable for the bench alignment of modular leg prosthesis systems. For the three-dimensional alignment of a modular leg prosthesis system, the knee joint, prosthetic foot and components are assembled according to the alignment recommendations. The socket is then secured in the device and attached to the knee joint.

1.3 Function

The PROS.A. Assembly was designed for the bench alignment of modular leg prosthesis systems (transfemoral and knee disarticulation prostheses). The knee joint is secured in the device according to the alignment values by using special adapter inserts, and is connected to the prosthetic foot and corresponding adapters (Figure 1–3). A socket clamping fixture is used for the three-dimensional positioning of the socket. The pivot point used to adjust the socket clamping fixture is aligned with the hip pivot point. This facilitates the reproducible alignment and fixation of the socket according to the patient's residual limb position (see Section 5 for particularities in regards to knee disarticulation prostheses).

The PROS.A. Assembly alignment apparatus facilitates the measurable, reproducible three-dimensional alignment of Ottobock modular leg prosthesis components without the use of a laser. All alignment recommendations are automatically taken into account if the basic settings are maintained. The finished modular leg prosthesis system can be measured, e.g. for subsequent fittings.

1.4 Connection Systems

Two versions of the socket clamping fixture of different diameters and with extensions as well as adapter inserts for Ottobock knee joints are included with the alignment apparatus (Figure 3 and 4). The 743Y594 socket clamping fixture is suitable for sockets with a minimum diameter of 145 mm and up to 200 mm; the 743Y594=1 socket clamping fixture is suitable for sockets with a minimum diameter of 90 mm and up to 145 mm.

1.5 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

All repair and service work must be performed by authorised personnel and under observation of the applicable legal regulations! Improper operation or repairs that were not completed by Ottobock may pose a risk to the user. The manufacturer is not liable for damage caused by improper use.

2 Safety Instructions

2.1 Explanation of Symbols

 WARNING	Warnings regarding possible risks of severe accident or injury.
 CAUTION	Warnings regarding possible risks of accident or injury.
 NOTICE	Warnings regarding possible technical damage.
 INFORMATION	Information regarding operation. Information for service personnel.

2.2 General Safety Instructions

INFORMATION

Please read the instructions for use first! Familiarise yourself with the handling and functionality of the product prior to use.

CAUTION

Risk of tripping / contact with protruding parts. The PROS.A. Assembly should not be installed free-standing but always with the back to a wall. Socket clamping fixtures always have to be stored on the adapter plate without the extension and with the clamping fixture pointing inward. The height support adjustment crank must be folded.

CAUTION

Risk of injury due to improperly clamped prosthesis components. Only prosthesis components that can be securely held in place and handled during the alignment of a transfemoral prosthesis may be clamped in the device.

⚠ CAUTION

Risk of injury as a result of improper use. Improper use is associated with the risk of incorrect alignment results and therefore the risk of injury for the patient. Only utilise the product in accordance with its intended use.

Compliance with the alignment instructions in Section 4.5 and in the instructions for use for the respective prosthesis joints / components is mandatory.

Installing the apparatus on a level surface is essential for precise alignment results and measurements.

NOTICE

Risk of damage due to severe vibrations. Exposing the PROS.A. Assembly to severe vibrations can cause structural distortion, leading to incorrect alignment results and measurements. The PROS.A. Assembly must not be exposed to severe vibrations.

NOTICE

Risk of damage to the hydraulics. The hydraulic adjustment of the socket adapter height support must not be turned past the stop.

NOTICE

Risk of damage to the socket clamping fixture. Sharp, pointed objects must be kept away from the inflatable components of the socket clamping fixture.

NOTICE

Risk of damage caused by incorrect packaging. Only use the original packing for shipping or transport.

INFORMATION

The socket adapter height support hydraulics are designed to be maintenance-free.

3 Delivery and Preparation for Use

3.1 Instructions for Safe Transportation

NOTICE

Risk of damage due to improper transportation. The PROS.A. Assembly alignment apparatus must not be transported without packaging. Severe vibrations during transportation must be avoided. Do not lift or transport the alignment apparatus by the knee joint adapter or socket adapter.

The PROS.A. Assembly is delivered in a wooden crate. Transport straps are used to secure the PROS.A. Assembly for transportation so that the risk of sliding and / or tipping is eliminated.

3.2 Unpacking and Checking

The PROS.A. Assembly is supplied fully assembled.

1. After unpacking, check the PROS.A. Assembly and accessories for completeness and possible damage incurred during transportation. Any transport damage must be reported to the carrier promptly.
2. Remove transportation safety devices, packaging materials and adhesives from the PROS.A. Assembly. Take the instructions for use and other accessories out of the package.

3.3 Place of Installation and Environmental Conditions

The ambient temperature should not significantly exceed 22 °C. The specified technical data (see Section 7) refer to an ambient temperature of 22 °C. The data may vary depending on the environmental conditions.

3.4 Installing the Alignment Apparatus

NOTICE

Risk of damage due to severe vibrations. Exposing the PROS.A. Assembly to severe vibrations can cause structural distortion, leading to incorrect alignment results and measurements. The PROS.A. Assembly must not be exposed to severe vibrations.

Installing the PROS.A. Assembly on a level surface is essential for precise alignment results and measurements. The height-adjustable feet (Figure 1, Number 11) are used for the **vertical alignment** of the apparatus. A laser (see Section 4.5) or level may be used to verify proper alignment.

4 Alignment of a Modular Transfemoral Prosthesis

INFORMATION

The PROS.A. Assembly is also suitable for the alignment of knee disarticulation prostheses. Please see the related information in Section 5.



Systematic prosthesis alignment is essential for correct prosthesis statics. Bench alignment guidelines established in clinical studies form the basis for proper alignment. Ottobock developed the alignment recommendations for modular transfemoral leg prosthesis systems based on these guidelines (see enclosed poster "Alignment Recommendations for Modular Transfemoral Leg Prosthesis Systems", article no. 646F219=*).


The design of the PROS.A. Assembly simplifies the implementation of these alignment recommendations. With the correct bench alignment of modular transfemoral leg prosthesis systems in the PROS.A. Assembly, the advantages of the knee joints are fully realised.

4.1 Foot Positioning

A recommended prosthetic foot is selected according to MOBIS® (selection based on the mobility grade and functional patient requirements).

1. Loosen the clamping screw below the footplate (Figure 5, Item A) and set the required alignment data – anterior placement, outward rotation of the foot and heel height at the foot adapter – as follows (Figure 5/6):
 - Select the anterior placement of the foot centreline in reference to the alignment reference line according to the information in the instructions for use for the corresponding foot (Figure 5, Item B). Please consult the instructions for use of the corresponding foot.
 - Customise the outward rotation of the foot. The default setting is approximately 5° (Figure 6). Tighten the clamping screw.
2. Set the effective heel height plus 5 mm (Figure 5, Item C).
3. Position the foot so the longitudinal axis is centred on the footplate. The centrelines of the heel and heel block have to be aligned. (Figure 7/8).

4.2 Knee Joint Positioning

⚠ CAUTION



Risk of pinching due to improper height support clamping. After each adjustment of the knee joint adapter height support, the clamping levers must be tightened to prevent the height support from moving down.

For the alignment of a modular transfemoral prosthesis in the PROS.A. Assembly, the position of the knee joint according to the Ottobock alignment recommendations is automatically specified by the knee bracket.

1. The outward rotation of 5° for the respective amputation side (left / right) is set using the knee bracket on the knee joint adapter height support. In order to do so, insert the pin of the knee bracket into the knee bracket adapter (Figure 9).

2. Insert the correct adapter inserts for the joint into the knee bracket. They are identified by one or two dots (Figure 10).
3. Clamp the modular knee joint in place. The outward rotation (Figure 11) as well as the frontal and sagittal position are precisely specified by the adapter insert and the knee bracket.
4. Determine the length of the tube adapter:
 - Push the height support up and secure it using the clamping levers.
 - Attach the tube clamp adapter to the tube adapter.
 - Loosen the adjustment screws on the tube adapter and attach the tube adapter to the pyramid adapter on the foot.
 - Move the height support down until the knee joint pyramid adapter engages in the tube clamp adapter (Figure 12). Read the knee-floor measurement from the scale (Figure 13).
 - Calculate the difference compared to the desired length and adjust the tube adapter as required.
5. Finally, mount the modular components using the adjustment screws.



4.3 Socket Positioning


4.3.1 Socket Clamping Fixture Functions

The PROS.A. Assembly is equipped with a socket quick-clamping fixture for socket positioning (Figure 15/16).

In order to customise the prosthesis to patient requirements, this fixture is used to position the socket against the alignment reference line according to the patient data as follows:

- **Flexion / extension:** ° setting on the pivot point in the socket clamping fixture (Figure 14, Item A). A clamping screw is used for fixation (Figure 14, Item B).
- **Adduction / abduction:** mm setting on the pivot point of the socket clamping fixture adapter (Figure 15, Item E). Basic setting (vertical position): **The lengthwise mark on the socket clamping fixture axle and centre mark on the socket clamping fixture adapter are aligned.** (unit°; 1 scale line 2°) The clamping lever is used for loosening / fixation in order to change the settings (Figure 15, Item B).
- **Anterior / posterior:** Anterior / posterior placement of the socket clamping fixture axle on the pivot point of the socket clamping fixture adapter (Figure 16, Item A). Basic setting: **2nd axle mark (recessed groove)** aligned flush with the socket clamping fixture adapter (Figure 15, Item C). The clamping lever is used for loosening / fixation (Figure 15, Item B).
- **Medial / lateral:** Medial or lateral positioning of the socket clamping fixture adapter on the profile tube (Figure 15, Item D). Basic setting (home position): **The socket clamping fixture adapter is centred over the M-L scale** (Figure 15, Item E). The clamping lever is used for loosening / fixation (Figure 15, Item F).
- The inflatable section of the socket clamping fixture (Figure 16, Item B), the bellows (Figure 16, Item C) and the centring unit with flexible star (Figure 16, Item D) are used to secure the socket.

4.3.2 Preparing the Socket

Prior to positioning / securing the socket, the socket centre must be marked with a proximal and distal point on the lateral side as follows:

- Distal point (Figure 17, Item A): Mark the socket centre on the end of the socket. Then draw a line through both points from the edge of the socket to the end.
- Proximal point (Figure 17, Item B): Mark the top edge of the trochanter major in the centre of the socket (hip joint pivot point).

4.3.3 Securing the Socket

INFORMATION


The socket is always positioned in the same manner, regardless of the socket material, technique, residual limb length and connection adapters.

Before positioning the socket, the hip joint pivot point (trochanter major mark) must be horizontally and vertically aligned with the socket clamping fixture pivot point and secured in this position.

The socket is centred by the centring unit with flexible star at the distal point (Figure 16, Item D) and the self-centring inflatable socket clamping fixture at the proximal point (Figure 16, Item B).

1. Select the correct socket clamping fixture according to the socket diameter (743Y594=1 for 90–145 mm, 743Y594 for 145–200 mm) as well as the right extension for the centring unit according to the socket length. Insert the socket clamping fixture into the adapter and secure it with the clamping lever (Figure 15, Item B).
2. Push the socket over the centring unit. The flexible star is used as a guide for centring at the distal end (Figure 18, Item A).

3. Press the socket against the spring-loaded centring unit until the marked point on the socket at the level of the trochanter major (Figure 19, Item B) is aligned with the socket clamping fixture pivot point (Figure 18/19, Item B). At the same time, ensure the socket is secured in the centred position at the distal end.
4. Close the screw valve on the bellows and inflate the inflatable section of the socket clamping fixture on the proximal end of the socket.
5. Check the pivot point position and correct if required.

4.3.4 Positioning the Socket

Now the socket clamping fixture adjustment options (see Section 4.3.1) are used to correctly position the socket in relation to the knee joint. The flexion / extension, abduction / adduction, A/P offset and rotation settings are based on the patient data and the specifications in the Ottobock alignment recommendations.

4.3.5 Connecting the Socket to the Knee Joint

1. In order to obtain assembly clearance, shift the socket clamping fixture in the proximal direction using the hydraulic socket adapter height support adjustment mechanism (Figure 19, Item C).
2. Install the rotation adapter / double adapter and lamination anchor or socket adapter + socket attachment block on the knee joint.
3. Distally lower the socket over the anchor arms or the socket attachment block using the hydraulic adjustment mechanism of the height support.
4. Adjust the “tuberosity-floor measurement” (Figure 1, Item 11) in accordance with the patient measurements using the hydraulic adjustment mechanism of the height support (Figure 20).
5. Connect the socket to the lamination anchor or socket attachment block.


4.4 Documenting the Prosthesis Alignment

The three-dimensional alignment of the modular leg prosthesis with the PROS.A. Assembly is measurable and reproducible.

The following alignment data can be read from the integrated scales (Figure I-IV) and used for subsequent fittings:

- **Scale I** (on the footplate): Foot position + outward rotation of the foot (Figure 21)
- **Scale II** (on the lateral frame): Knee-floor measurement (Figure 13); trochanter major floor measurement (Figure 22)
- **Scale III** (on the socket clamping fixture adapter pivot point): Adduction / abduction + A-P position of the socket (Figure 23)
- **Scale IV** (on the socket clamping fixture pivot point): Socket flexion / extension (Figure 24)


4.5 Using a Laser

The alignment of a modular leg prosthesis system in reference to the alignment reference line using the PROS.A. Assembly can be verified using an external laser (e.g. 743L5 crosshair laser or 743L20=230 LaserLine). The alignment apparatus and laser must be aligned as follows:

1. Align the PROS.A. Assembly with crosshairs drawn on the floor (or existing tile joints at right angles).
2. Use the height-adjustable feet to position the alignment apparatus so it is vertical (Figure 1, Item 12). Verify the vertical A-P and frontal alignment using a laser plumb if applicable. Use the lateral frame of the PROS.A. Assembly as a reference line (Figure 25).

When verifying prosthesis alignment using a laser, ensure that the laser plumb line is **precisely vertically aligned** with the respective arrows on the frame (Figure 26). The arrows identify the alignment reference line (Figure 26, Item A; Figure 27) and the frontal plumb line (Figure 26, Item B; Figure 28).


4.6 Static and Dynamic Alignment Optimisation


After the bench alignment of the prosthesis components on the PROS.A. Assembly is complete according to the alignment recommendations, the static alignment optimisation of the transfemoral prosthesis on the L.A.S.A.R. Posture (Figure 29) follows in accordance with the “Alignment Recommendations for Modular Transfemoral Leg Prosthesis Systems”, article no. 646F219=* by adjusting the plantar flexion position of the foot under load.

The systematic approach described above correctly aligns the prosthesis using a specific, reproducible process.

The prosthetist verifies the settings based on the movement pattern during dynamic trial walking.

5 Alignment of Knee Disarticulation Prostheses

The PROS.A. Assembly is also suitable for the alignment of knee disarticulation prostheses. The following particularities must be taken into account:

Due to the limited socket length, the socket clamping fixture pivot point cannot be aligned with the hip pivot point (trochanter major).

The difference between the trochanter major height and socket clamping fixture pivot point height must be compensated through the anterior placement of the socket clamping fixture in relation to the alignment reference line as the flexion position of the socket increases.

Procedure: After securing the socket in the alignment apparatus, the height of the quick-clamping fixture pivot point is read from scale II (Figure 22).

The values for the required anterior placement of the socket clamping fixture at various socket flexion angles can be read from the table to the left according to the calculated difference between the height of the trochanter major and the height of the socket clamping fixture pivot point.

For anterior placement, the axle of the socket clamping fixture is loosened at the socket clamping fixture adapter pivot point and moved in the anterior direction (see Section 4.3.1; “Anterior / posterior”).

INFORMATION The maximum displacement of the socket clamping fixture axle is 30 mm anterior and 10 mm posterior. The distance from groove to groove on the axle is 5 mm.

Socket flexion (in °)				Distance between the trochanter major and the socket adapter pivot point (in mm)
15	12	7	3	
7	5	3	1	25
13	10	6	3	50
20	16	9	4	75
27	21	12	5	100
--	26	15	7	125
Compensation by anterior placement of the socket clamping fixture axle (in mm)				

6 Maintenance and Cleaning

6.1 Cleaning

NOTICE

Risk of damage due to improper cleaning. Do not use any aggressive cleaners, solvents, fuels, acids or hard brushes etc. Sharp, pointed objects must be kept away from the inflatable components of the socket clamping fixture.

Regular general cleaning of the PROS.A. Assembly helps maintain the trouble-free operation of the height-adjustable supports for the knee and socket adapter. It also extends the life of the PROS.A. Assembly.

Use a damp cloth and mild cleaning solution or, if required, a moist plastic brush to remove dirt.

7 Technical Data

Dimensions and Weights	
Dimensions (WxHxD)	700 mm × 2000 mm × 680 mm
Gross weight (without prosthesis)	45 kg
Maximum prosthesis height	1200 mm
Environmental Conditions	
Allowable temperature range – storage	- 10 °C ... + 60 °C
Allowable temperature range – operation	+ 5 °C ... + 40 °C
Allowable humidity during operation and storage	70% relative humidity, non-condensing

8 Disposal Information

If the PROS.A. Assembly alignment apparatus is no longer used, disposal by a local disposal company in accordance with the applicable environmental protection laws in the respective country is required.

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Otto Bock HealthCare GmbH