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SIL

Silicone

is durable and easy to clean. The material provides high stability and good adhesion. Silicone is suitable for fittings on residual limbs with good soft tissue coverage.

This material is used for both transtibial and transfemoral prostheses. We recommend silicone liners for patients wearing a transtibial prosthesis, with a low to moderate mobility grade. Patients with a transfemoral prosthesis can wear a SIL liner regardless of their mobility grade.

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
TPE

Copolymer

is highly elastic. It contains skin-friendly mineral white oil, making it especially suitable for residual limbs with dry skin. We recommend copolymer liners for patients with a low to medium mobility grade, insofar as they wear a transtibial prosthesis. This liner material is rarely used for transfemoral amputees.

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
PUR

Polyurethane

ensures even pressure distribution and offers a precise, comfortable fit. Optimum dampening of impacts also makes the material suitable for sensitive, bony and/or scarred residual limbs. We recommend polyurethane for patients with a low to especially high mobility grade, regardless of the amputation level.

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
TT fitting suggestions based on different residual limb types

Residual limb shapes versus mobility grade

Cylindrical residual limb shape with very tender soft tissue coverage and overhang

	Skeo		Skeo		Skeo		Skeo	
Pin	Caleo		Caleo		Caleo		Caleo	
Valve	6Y41		6Y41		6Y41		6Y41	
DVS	6Y92/95=C		6Y92/95=C		6Y92/95=C		6Y92/95=C	
Harmony	6Y523/Unique		6Y523/Unique		6Y523/Unique		6Y523/Unique	
	Skeo		Skeo		Skeo		Skeo	
	Caleo		Caleo		Caleo		Caleo	
	6Y94		6Y94		6Y94		6Y94	
	Unique 맞춤형		Unique 맞춤형		Unique 맞춤형		Unique 맞춤형	

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TT fitting suggestions based on different residual limb types

Residual limb shapes versus mobility grade

Conical residual limb shape with tender soft tissue coverage and overhang


Pin 	Skeo Caleo Uneo	Skeo Caleo Uneo	Skeo Caleo Uneo	Skeo Caleo Uneo
Valve 	Skeo 6Y41 Caleo 6Y92/95=C Uneo	Skeo 6Y41 Caleo 6Y92/95=C Uneo 6Y512/Unique	Skeo 6Y41 Caleo 6Y92/95=C Uneo 6Y512/Unique	Skeo 6Y41 Caleo 6Y92/95=C Uneo 6Y512/Unique
DVS 	Skeo Caleo Uneo	Skeo Caleo 6Y94 Uneo	Skeo Caleo 6Y94 Uneo	Skeo Caleo 6Y94 Uneo
Harmony 	Skeo Caleo Uneo	Skeo Caleo Uneo 6Y512/Unique	Skeo Caleo Uneo 6Y512/Unique	Skeo Caleo Uneo 6Y512/Unique

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TT fitting suggestions based on different residual limb types

Residual limb shapes versus mobility grade

Conical residual limb shape with good soft tissue coverage

Pin	Skeo 6Y43 Caleo 6Y95=L Uneo	Skeo 6Y43 Caleo 6Y95=L Uneo	Skeo 6Y43 Caleo Uneo	Skeo Caleo Uneo
Valve	Skeo Caleo 6Y95=C Uneo	Skeo 6Y41 Caleo 6Y95=C Uneo	Skeo 6Y41 Caleo 6Y95=C Uneo 6Y512/Unique	Skeo 6Y41 Caleo 6Y95=C Uneo 6Y512/Unique
DVS	Skeo Caleo Uneo	Skeo Caleo 6Y94 Uneo	Skeo Caleo 6Y94 Uneo	Skeo Caleo 6Y94 Uneo
Harmony	Skeo Caleo Uneo	Skeo Caleo Uneo 6Y512/Unique	Skeo Caleo Uneo 6Y512/Unique	Skeo Caleo Uneo 6Y512/Unique

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TT fitting suggestions based on different residual limb types

Residual limb shapes versus mobility grade

Bony residual limb shape

Pin		Skeo Caleo 6Y95=L Uneo	Skeo Caleo 6Y95=L Uneo	Skeo Caleo Uneo	Skeo Caleo Uneo
Valve		Skeo Caleo 6Y95=C Uneo	Skeo Caleo 6Y92/95=C Uneo 6Y512/Unique	Skeo Caleo 6Y92/95=C Uneo 6Y512/Unique	Skeo Caleo 6Y92/95=C Uneo 6Y512/Unique
DVS		Skeo Caleo Uneo	Skeo Caleo 6Y94 Uneo	Skeo Caleo 6Y94 Uneo	Skeo Caleo 6Y94 Uneo
Harmony		Skeo Caleo Uneo	Skeo Caleo Uneo 6Y512/Unique	Skeo Caleo Uneo 6Y512/Unique	Skeo Caleo Uneo 6Y512/Unique

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TT fitting suggestions based on different residual limb types

Residual limb shapes versus mobility grade

Sensitive residual limbs with retracted scars

Pin	Skeo	Skeo	Skeo	Skeo
	Caleo	Caleo	Caleo	Caleo
	Uneo	Uneo	Uneo	Uneo
Valve	Skeo	Skeo	Skeo	Skeo
	Caleo 6Y92/95=C	Caleo 6Y92/95=C	Caleo 6Y92/95=C	Caleo 6Y92/95=C
	Uneo	Uneo 6Y512/Unique	Uneo 6Y512/Unique	Uneo 6Y512/Unique
DVS	Skeo	Skeo	Skeo	Skeo
	Caleo	Caleo 6Y94	Caleo 6Y94	Caleo 6Y94
	Uneo	Uneo	Uneo	Uneo
Harmony	Skeo	Skeo	Skeo	Skeo
	Caleo	Caleo	Caleo	Caleo
	Uneo	Uneo 6Y512/Unique	Uneo 6Y512/Unique	Uneo 6Y512/Unique

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TF fitting suggestions based on different residual limb types

Cylindrical residual limb shape with very tender soft tissue coverage and overhang

Residual limb shapes versus mobility grade

	Mobility grade 1		Mobility grade 2		Mobility grade 3		Mobility grade 4	
KISS	Skeo 6Y80/85		Skeo 6Y80/85		Skeo 6Y80/85		Skeo	
	Caleo		Caleo		Caleo		Caleo	
	Uneo		Uneo		Uneo		Uneo	
Valve + ProSeal	Skeo		Skeo		Skeo 6Y110/81		Skeo 6Y110/81	
	Caleo		Caleo		Caleo		Caleo	
	Uneo		Uneo		Uneo		Uneo	
Harmony	Skeo		Skeo		Skeo 6Y110/81		Skeo 6Y110/81	
	Caleo		Caleo		Caleo		Caleo	
	Uneo		Uneo		Uneo		Uneo	

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

Tested according to JIS Z 2801 with *Staphylococcus aureus*

Germ count

■ Bacterial growth on untreated polymer

■ Skinguard inhibits bacterial growth

○ Skinguard inhibits the growth of bacteria on polymer surfaces.

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